

兔抗 TAFA4 多克隆抗体

中文名称: 兔抗 TAFA4 多克隆抗体

英文名称: Anti-TAFA4 rabbit polyclonal antibody

别 名: TAFA chemokine like family member 4; TAFA-4; FAM19

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗原: TAFA4

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

This gene is a member of the TAFA family which is com posed of five highly homologous genes that encode sma II secreted proteins. These proteins contain conserved cys teine residues at fixed positions, and are distantly related to MIP-1alpha, a member of the CC-chemokine family. The TAFA proteins are predominantly expressed in specific regions of the brain, and are postulated to function as brain-specific chemokines or neurokines, that act as regula tors of immune and nervous cells. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Nov 2011]



Applications:	ELISA, WB, IHC
Name of antibody:	TAFA4
Immunogen:	Fusion protein of human TAFA4
Full name:	TAFA chemokine like family member 4
Synonyms:	TAFA-4; FAM19A4
SwissProt:	Q96LR4
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	16 kDa
WB Positive control:	Rat kidney tissue lysate
WB Recommended dilution:	500-2000





