

兔抗 LRP3 多克隆抗体

中文名称: 兔抗 LRP3 多克隆抗体

英文名称: Anti-LRP3 rabbit polyclonal antibody

相关类别: 一抗

储 存: 冷冻 (-20℃)

宿 主: Rabbit

抗 原: LRP3

反应种属: Human

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Low density lipoprotein receptor-related protein 3 (LRP-3) is a protein that in humans is encoded by theLRP3 gene. Probable receptor, which may be involved in the internalization of lipophilic molecules and/or signal transduction. Its precise role is however unclear, since it does not bind to very low density lipoprotein (VLDL) or to LRPAP1 in vitro.
Applications:	ELISA, IHC
Name of antibody:	LRP3
Immunogen:	Synthetic peptide of human LRP3
Full name:	low density lipoprotein receptor-related protein 3
SwissProt:	O75074
ELISA Recommended dilutio	2000-10000

n:	
IHC positive control:	Human brain
IHC Recommend dilution:	100-300

