

小鼠抗 ETS2 单克隆抗体

中文名称: 小鼠抗 ETS2 单克隆抗体

英文名称: Anti-ETS2 mouse monoclonal antibody

别名: ETS2IT1

储存: 冷冻 (-20℃) 避光

抗原: ETS2

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	ETS transcriptions factors, such as ETS2, regulate numerous genes and are involved in stem cell development, cell senescence and death, and tumorigenesis. The conserved ETS domain within these proteins is a winged helix-turn-helix DNA-binding domain that recognizes the core consensus DNA sequence GGAA/T of target genes (Dwyer et al., 2007 [PubMed 17986575]).[supplied by OMIM].
Applications:	WB, IHC, IF
Name of antibody:	ETS2
Immunogen:	Fusion protein of human ETS2

Full name:	v-ets avian erythroblastosis virus E26 oncogene homolog 2 (ETS2), transcript variant 1
Synonyms:	ETS2IT1
SwissProt:	P15036
IHC positive control:	human bladder tissue and carcinoma of human bladder tissue; human endometrium tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	53 kDa
WB Positive control:	HepG2 and Hela cell lysates
WB Recommended dilution:	200-1000