

## 兔抗 LRP10 多克隆抗体

中文名称: 兔抗 LRP10 多克隆抗体

英文名称: Anti-LRP10 rabbit polyclonal antibody

别名: LRP9; MST087; MSTP087

抗原: LRP10

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

相关类别: 一抗

反应种属: Human Mouse

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	Probable receptor, which is involved in the internalization of lipophilic molecules and/or signal transduction. May be involved in the uptake of lipoprotein APOE in liver By similarity.
<b>Applications:</b>	IHC
<b>Name of antibody:</b>	LRP10
<b>Immunogen:</b>	Synthesized peptide derived from internal of human LRP10.
<b>Full name:</b>	low density lipoprotein receptor-related protein 10
<b>Synonyms:</b>	LRP9; MST087; MSTP087
<b>SwissProt:</b>	Q7Z4F1

<b>IHC positive control:</b>	Human brain tissue
<b>IHC Recommendation dilution:</b>	50-100

