

兔抗 LRP11 多克隆抗体

中文名称：兔抗 LRP11 多克隆抗体

英文名称：Anti-LRP11 rabbit polyclonal antibody

别名：MANSC3; bA350J20.3

抗原：LRP11

储存：冷冻（-20℃）避光

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Belongs to the LDLR family. 2 isoforms of the human protein are produced by alternative splicing.
Applications:	WB
Name of antibody:	LRP11
Immunogen:	Synthesized peptide derived from internal of human LRP11.
Full name:	low density lipoprotein receptor-related protein 11
Synonyms :	MANSC3; bA350J20.3
SwissProt:	Q86VZ4
WB Predicted band size	53 kDa

e:	
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

