

兔抗 TAX1BP1 多克隆抗体

中文名称：兔抗 TAX1BP1 多克隆抗体

英文名称： Anti-TAX1BP1 rabbit polyclonal antibody

别名： T6BP; TXBP151; CALCOCO3

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TAX1BP1

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a HTLV-1 tax1 binding protein. The encoded protein interacts with TNFAIP3, and inhibits TNF-induced apoptosis by mediating the TNF AIP3 anti-apoptotic activity. Degradation of this protein by caspase-3-like family proteins is associated with apoptosis induced by TNF. This protein may also have a role in the inhibition of inflammatory signaling pathways. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC

Name of antibody:	TAX1BP1
Immunogen:	Fusion protein of human TAX1BP1
Full name:	Tax1 (human T-cell leukemia virus type I) binding protein 1
Synonyms :	T6BP; TXBP151; CALCOCO3
SwissProt:	Q86VP1
ELISA Recommended dilution:	3000-10000
IHC positive control:	Human colon cancer and human esophagus cancer
IHC Recommend dilution:	100-300

