

兔抗 TAZ 多克隆抗体

- 中文名称: 兔抗 TAZ 多克隆抗体
- 英文名称: Anti-TAZ rabbit polyclonal antibody
- 别 名: EFE; BTHS; EFE2; G4.5; Taz1; CMD3A; LVNCX
- 相关类别: 一抗
- 储 存: 冷冻 (-20℃) 避光
- 宿 主: Rabbit
- 抗 原: TAZ
- 反应种属: Human
- 标记物: Unconjugate
- 克隆类型: Unconjugate

技术规格

	This gene encodes a protein that is expressed at high lev
	els in cardiac and skeletal muscle. Mutations in this gene
	have been associated with a number of clinical disorders
	including Barth syndrome, dilated cardiomyopathy (DCM),
	hypertrophic DCM, endocardial fibroelastosis, and left ven
Background:	tricular noncompaction (LVNC). Multiple transcript variants
	encoding different isoforms have been described. A long
	form and a short form of each of these isoforms is prod
	uced; the short form lacks a hydrophobic leader sequenc
	e and may exist as a cytoplasmic protein rather than bei
	ng membrane-bound. Other alternatively spliced transcript



	s have been described but the full-length nature of all th ese transcripts is not known.
Applications:	WB
Name of antibody:	TAZ
Immunogen:	Synthesized peptide derived from Internal of human TAZ.
Full name:	tafazzin
Synonyms :	EFE; BTHS; EFE2; G4.5; Taz1; CMD3A; LVNCX
SwissProt:	Q16635
WB Predicted band size:	33 kDa
WB Positive control:	HepG2 cells lysate
WB Recommended dilution:	500-3000

