

小鼠抗 FCMR 单克隆抗体

中文名称: 小鼠抗 FCMR 单克隆抗体

英文名称: Anti-FCMR mouse monoclonal antibody

别名: TOSO; FAIM3

抗原: FCMR

储存: 冷冻 (-20°C) 避光

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	May play a role in the immune system processes. Protects cells from FAS-, TNF alpha- and FADD-induced apoptosis without increasing expression of the inhibitors of apoptosis BCL2 and BCLXL. Seems to activate an inhibitory pathway that prevents CASP8 activation following FAS stimulation, rather than blocking apoptotic signals downstream. May inhibit FAS-induced apoptosis by preventing CASP8 processing through CFLAR up-regulation.
Applications:	IHC
Name of antibody:	FCMR
Immunogen:	Fusion protein of human FCMR

Full name:	Fc fragment of IgM receptor (FCMR), transcript variant 1
Synonyms:	TOSO; FAIM3
SwissProt:	O60667
IHC positive control:	carcinoma of human thyroid tissue
IHC Recommend dilution:	30-150