

## 兔抗 TBC1D2 多克隆抗体

中文名称：兔抗 TBC1D2 多克隆抗体

英文名称： Anti-TBC1D2 rabbit polyclonal antibody

别名： TBC1 domain family member 2; PARIS1; PARIS-1; TBC1D2A

相关类别： 一抗 储存 冷冻 (-20℃)

宿主： Rabbit

抗原： TBC1D2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Acts as GTPase-activating protein for RAB7A. Signal effector acting as a linker between RAC1 and RAB7A, leading to RAB7A inactivation and subsequent inhibition of cadherin degradation and reduced cell-cell adhesion.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	TBC1D2
<b>Immunogen:</b>	Fusion protein of human TBC1D2
<b>Full name:</b>	TBC1 domain family member 2
<b>Synonyms:</b>	PARIS1; PARIS-1; TBC1D2A
<b>SwissProt:</b>	Q9BYX2
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	30-150

