

小鼠抗 FERMT3 单克隆抗体

中文名称： 小鼠抗 FERMT3 单克隆抗体

英文名称： Anti-FERMT3 mouse monoclonal antibody

别名： fermitin family member 3; URP2; KIND3; MIG-2; MIG2B; URP2SF; UNC112C

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： FERMT3

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	Kindlins are a small family of proteins that mediate protein-protein interactions involved in integrin activation and thereby have a role in cell adhesion, migration, differentiation, and proliferation. The protein encoded by this gene has a key role in the regulation of hemostasis and thrombosis. This protein may also help maintain the membrane skeleton of erythrocytes. Mutations in this gene cause the autosomal recessive leukocyte adhesion deficiency syndrome-III (LAD-III). Alternative splicing results in multiple transcript variants encoding distinct isoforms.
Applications:	WB, IHC

Name of antibody:	FERMT3
Immunogen:	Fusion protein of human FERMT3
Full name:	fermitin family member 3
Synonyms:	URP2; KIND3; MIG-2; MIG2B; URP2SF; UNC112C
SwissProt:	Q86UX7
IHC positive control:	Human tonsil tissue
IHC Recommend dilution:	50-150
WB Predicted band size:	76 KD
WB Positive control:	Jurkat cell lysate
WB Recommended dilution:	200-1000