

小鼠抗 FGG 单克隆抗体

中文名称: 小鼠抗 FGG 单克隆抗体

英文名称: Anti-FGG mouse monoclonal antibody

别名: fibrinogen, gamma chain; fibrinogen, gamma polypeptide; fibrinogen gamma chain

抗原: FGG

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:

The protein encoded by this gene is the gamma component of fibrinogen, a blood-borne glycoprotein comprised of three pairs of nonidentical polypeptide chains. Following vascular injury, fibrinogen is cleaved by thrombin to form fibrin which is the most abundant component of blood clots. In addition, various cleavage products of fibrinogen and fibrin regulate cell adhesion and spreading, display vasoconstrictor and chemotactic activities, and are mitogens for several cell types. Mutations in this gene lead to several disorders, including dysfibrinogenemia, hypofibrinogenemia and thrombophilia. Alternative splicing results in two transcript variants encoding different isoforms. [provided by RefSeq, Jul

Applications:	WB, IF
Name of antibody:	FGG
Immunogen:	Fusion protein of human FGG
Full name:	fibrinogen gamma chain (FGG), transcript variant gamma-A
Synonyms:	fibrinogen, gamma chain; fibrinogen, gamma polypeptide; fibrinogen gamma chain
SwissProt:	P02679
WB Predicted band size:	52 kDa
WB Positive control:	HepG2 cell lysate
WB Recommended dilution:	500-2000