

兔抗 LRRC28 多克隆抗体

中文名称: 兔抗 LRRC28 多克隆抗体

英文名称: Anti-LRRC28 rabbit polyclonal antibody

别名: leucine rich repeat containing 28

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: LRRC28

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The function of the LRRC28 protein is not known.
Applications:	ELISA, WB, IHC
Name of antibody:	LRRC28
Immunogen:	Fusion protein of human LRRC28
Full name:	leucine rich repeat containing 28
SwissProt:	Q86X40
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human brain
IHC Recommend dilution:	50-200

WB Predicted band size:	42 kDa
WB Positive control:	Jurkat and HepG2 cell lysates
WB Recommended dilution:	500-2000





