

## 兔抗 LRRC41 多克隆抗体

中文名称：兔抗 LRRC41 多克隆抗体

英文名称：Anti-LRRC41 rabbit polyclonal antibody

别名：MUF1; PP7759

抗原：LRRC41

储存：冷冻（-20℃）避光

宿主：Rabbit

相关类别：一抗

反应种属：Human Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

|                          |  |
|--------------------------|--|
| <b>Background:</b>       | Probable substrate recognition component of an ECS (Elongin BC-CUL2/5-SOCS-box protein) E3 ubiquitin ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins. |
| <b>Applications:</b>     | WB   |
| <b>Name of antibody:</b> | LRRC41   |
| <b>Immunogen:</b>        | Synthesized peptide derived from internal of human LRRC41.   |
| <b>Full name:</b>        | leucine rich repeat containing 41  |
| <b>Synonyms :</b>        | MUF1; PP7759   |

|                                 |   |
|---------------------------------|---|
| <b>SwissProt:</b>               | Q15345                                  |
| <b>WB Predicted band size:</b>  | 89 kDa                                  |
| <b>WB Positive control:</b>     | Rat kidney and rat heart tissue lysates |
| <b>WB Recommended dilution:</b> | 500-3000                                |

