

兔抗 LRRC45 多克隆抗体

中文名称: 兔抗 LRRC45 多克隆抗体

英文名称: Anti-LRRC45 rabbit polyclonal antibody

别名: leucine rich repeat containing 45

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: LRRC45

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

| | |
|------------------------------------|--|
| Background: | Component of the proteinaceous fiber-like linker between two centrioles, required for centrosome cohesion. |
| Applications: | ELISA, WB, IHC |
| Name of antibody: | LRRC45 |
| Immunogen: | Fusion protein of human LRRC45 |
| Full name: | leucine rich repeat containing 45 |
| SwissProt: | Q96CN5 |
| ELISA Recommended dilution: | 5000-10000 |

| | |
|---------------------------------|--|
| IHC positive control: | Human esophagus cancer and Human colorectal cancer |
| IHC Recommend dilution: | 100-300 |
| WB Predicted band size: | 76 kDa |
| WB Positive control: | HepG2 cell lysate |
| WB Recommended dilution: | 500-2000 |





