

兔抗 LRRC49 多克隆抗体

中文名称: 兔抗 LRRC49 多克隆抗体

英文名称: Anti-LRRC49 rabbit polyclonal antibody

别名: leucine rich repeat containing 49; PGs4

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: LRRC49

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	LRRC49
Immunogen:	Synthetic peptide of human LRRC49
Full name:	leucine rich repeat containing 49
Synonyms:	PGs4
SwissProt:	Q8IUZ0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil

IHC Recommend dilution:

25-50

