

兔抗 TBPL1 多克隆抗体

中文名称：兔抗 TBPL1 多克隆抗体

英文名称： Anti-TBPL1 rabbit polyclonal antibody

别名： TATA-box binding protein like 1; TLF; TLP; STUD; TRF2; MGC:8389; MGC:9620

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TBPL1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a member of the TATA box-binding protein family. TATA box-binding proteins play a critical role in transcription by RNA polymerase II as components of the transcription factor IID (TFIID) complex. The encoded protein does not bind to the TATA box and initiates transcription from TATA-less promoters. This gene plays a critical role in spermatogenesis, and single nucleotide polymorphisms in this gene may be associated with male infertility. Alternatively spliced transcript variants have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 3.

Applications:	ELISA, IHC
Name of antibody:	TBPL1
Immunogen:	Fusion protein of human TBPL1
Full name:	TATA-box binding protein like 1
Synonyms:	TLF; TLP; STUD; TRF2; MGC:8389; MGC:9620
SwissProt:	P62380
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human lung cancer
IHC Recommend dilution:	40-200



