

小鼠抗 FLOT1 单克隆抗体

中文名称：小鼠抗 FLOT1 单克隆抗体

英文名称：Anti-FLOT1 mouse monoclonal antibody

别名：flotillin 1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：FLOT1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：mouse monoclonal

技术规格

Background:	This gene encodes an protein that localizes to the caveolae, which are small domains on the inner cell membranes. This protein plays a role in vesicle trafficking and cell morphology. Alternative splicing results in multiple transcript variants.
Applications:	WB
Name of antibody:	FLOT1
Immunogen:	Fusion protein of human FLOT1
Full name:	flotillin 1
SwissProt:	O75955
WB Predicted band size:	47 kDa

WB Positive control:	Mouse brain tissue lysate
WB Recommended dilution:	1000-5000

