

兔抗 LTBP4 多克隆抗体

中文名称：兔抗 LTBP4 多克隆抗体

英文名称：Anti-LTBP4 rabbit polyclonal antibody

抗 原：LTBP4

级 别：ARCL1C; LTBP-4; LTBP4L; LTBP4S

储 存：冷冻（-20℃）

宿 主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

The protein encoded by this gene binds transforming growth factor beta (TGFB) as it is secreted and targeted to the extracellular matrix. TGFB is biologically latent after secretion and insertion into the extracellular matrix, and sheds TGFB and other proteins upon activation. Defects in this gene may be a cause of cutis laxa and severe pulmonary, gastrointestinal, and urinary abnormalities. Three transcript variants enco

	ding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	LTBP4
Immunogen:	Synthetic peptide of human LTBP4
Full name:	latent transforming growth factor beta binding protein 4
Synonyms:	ARCL1C; LTBP-4; LTBP4L; LTBP4S
SwissProt:	Q8N2S1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and human breast cancer
IHC Recommend dilution:	25-100



