

小鼠抗 FOXO1(C-term)单克隆抗体

中文名称： 小鼠抗 FOXO1(C-term)单克隆抗体

英文名称： Anti-FOXO1(C-term) mouse monoclonal antibody

别名： forkhead box O1; FKH1; FKHR; FOXO1A

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： FOXO1(C-term)

反应种属： Human

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:	This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in myogenic growth and differentiation. Translocation of this gene with PAX3 has been associated with alveolar rhabdomyosarcoma.
Applications:	WB
Name of antibody:	FOXO1(C-term)
Immunogen:	Fusion protein of human FOXO1
Full name:	forkhead box O1
Synonyms:	FKH1; FKHR; FOXO1A

SwissProt:	Q12778
WB Predicted band size:	70 kDa
WB Positive control:	A431 and Hela cell lysates
WB Recommended dilution:	1000-5000

