

## 小鼠抗 FOXO1 单克隆抗体

中文名称： 小鼠抗 FOXO1 单克隆抗体

英文名称： Anti-FOXO1 mouse monoclonal antibody

别名： forkhead box O1; FKH1; FKHR; FOXO1A

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： FOXO1

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in myogenic growth and differentiation. Translocation of this gene with PAX3 has been associated with alveolar rhabdomyosarcoma.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	FOXO1
<b>Immunogen:</b>	Fusion protein of human FOXO1
<b>Full name:</b>	forkhead box O1
<b>Synonyms:</b>	FKH1; FKHR; FOXO1A
<b>SwissProt:</b>	Q12778
<b>WB Predicted band size:</b>	70 KD

<b>WB Positive control:</b>	Hela, 293T, HepG2 cell lysates
<b>WB Recommended dilution:</b>	500-2000

