

小鼠抗 FOXO4 单克隆抗体

中文名称： 小鼠抗 FOXO4 单克隆抗体

英文名称： Anti-FOXO4 mouse monoclonal antibody

别名： AFX; AFX1; MLLT7

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Mouse

抗原： FOXO4

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	This gene encodes a member of the O class of winged helix/forkhead transcription factor family. Proteins encoded by this class are regulated by factors involved in growth and differentiation indicating they play a role in these processes. A translocation involving this gene on chromosome X and the homolog of the Drosophila trithorax gene, encoding a DNA binding protein, located on chromosome 11 is associated with leukemia. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2010].
Applications:	WB, IHC

Name of antibody:	FOXO4
Immunogen:	Fusion protein of human FOXO4
Full name:	forkhead box O4 (FOXO4), transcript variant 1
Synonyms:	AFX; AFX1; MLLT7
SwissProt:	P98177
IHC positive control:	carcinoma of human pancreas tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	54 kDa
WB Positive control:	SVT2 cell lysate
WB Recommended dilution:	200-1000