

小鼠抗 FOXP1 单克隆抗体

中文名称: 小鼠抗 FOXP1 单克隆抗体

英文名称: Anti-FOXP1 mouse monoclonal antibody

别名: 12CC4; hFKH1B; HSPC215; MFH; QRF1

抗原: FOXP1

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:

This gene belongs to subfamily P of the forkhead box (FOX) transcription factor family. Forkhead box transcription factors play important roles in the regulation of tissue- and cell type-specific gene transcription during both development and adulthood. Forkhead box P1 protein contains both DNA-binding- and protein-protein binding-domains. This gene may act as a tumor suppressor as it is lost in several tumor types and maps to a chromosomal region (3p14.1) reported to contain a tumor suppressor gene(s). Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul

Applications:	WB, IHC, IF
Name of antibody:	FOXP1
Immunogen:	Fusion protein of human FOXP1
Full name:	forkhead box P1 (FOXP1), transcript variant 1
Synonyms:	12CC4; hFKH1B; HSPC215; MFH; QRF1
SwissProt:	Q9H334
IHC positive control:	human pancreas tissue and human lymphoma tissue; adenocarcinoma of human colon tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	75 kDa
WB Positive control:	Hela, A549, COS7, Jurkat, MCF-7 cell lysates
WB Recommended dilution:	200-1000