

小鼠抗 FOXP3 单克隆抗体

中文名称： 小鼠抗 FOXP3 单克隆抗体

英文名称： Anti-FOXP3 mouse monoclonal antibody

别名： forkhead box P3; JM2; AIID; IPEX; PIDX; XPID; DIETER

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： FOXP3

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	The protein encoded by this gene is a member of the forkhead/winged-helix family of transcriptional regulators. Defects in this gene are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX), also known as X-linked autoimmunity-immunodeficiency syndrome. Alternatively spliced transcript variants encoding different isoforms have been identified.
Applications:	ELISA, WB
Name of antibody:	FOXP3
Immunogen:	Fusion protein of human FOXP3
Full name:	forkhead box P3

Synonyms:	JM2; AIID; IPEX; PIDX; XPID; DIETER
SwissProt:	Q9BZS1
WB Predicted band size:	47 kDa
WB Positive control:	Mouse thymus tissue lysate
WB Recommended dilution:	500-2000

