

小鼠抗 FTCD 单克隆抗体

中文名称: 小鼠抗 FTCD 单克隆抗体

英文名称: Anti-FTCD mouse monoclonal antibody

别名: LCHC1

储存: 冷冻 (-20°C) 避光

抗原: FTCD

宿主: Mouse

反应种属: Human , Monkey

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	The protein encoded by this gene is a bifunctional enzyme that channels 1-carbon units from formiminoglutamate, a metabolite of the histidine degradation pathway, to the folate pool. Mutations in this gene are associated with glutamate formiminotransferase deficiency. Alternatively spliced transcript variants have been found for this gene.[provided by RefSeq, Dec 2009].
Applications:	WB, FC
Name of antibody:	FTCD
Immunogen:	Fusion protein of human FTCD
Full name:	formimidoyltransferase cyclodeaminase (FTCD), transcript variant

	t B
Synonyms:	LCHC1
SwissProt:	O95954
WB Predicted band size:	59 kDa
WB Positive control:	HepG2, COS7 cell lysates
WB Recommended dilution:	500-2000