

兔抗 TBX10 多克隆抗体

中文名称：兔抗 TBX10 多克隆抗体

英文名称： Anti-TBX10 rabbit ployclonal antibody

别名： TBX7; TBX13

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TBX10

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the T-box family of transcription factors. These transcription factors share a DNA-binding domain called the T-box, and play a role in several developmental processes including early embryonic cell fate and organogenesis. The encoded protein is a member of the T-box 1 subfamily. Mutations in this gene are thought to be a cause of isolated cleft lip with or without cleft palate.
Applications:	ELISA, IHC
Name of antibody:	TBX10

Immunogen:	Fusion protein of human TBX10
Full name:	T-box 10
Synonyms :	TBX7; TBX13
SwissProt:	O75333
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	25-100

