

小鼠抗 FUBP3 单克隆抗体

中文名称：小鼠抗 FUBP3 单克隆抗体

英文名称：Anti-FUBP3 mouse monoclonal antibody

别名：far upstream element binding protein 3; FBP3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：FUBP3

反应种属：Human

标记物：Unconjugate

克隆类型：mouse monoclonal

技术规格

Background:	May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression.
Applications:	WB
Name of antibody:	FUBP3
Immunogen:	Fusion protein of human FUBP3
Full name:	far upstream element binding protein 3
Synonyms:	FBP3
SwissProt:	Q96I24
WB Predicted band size:	62 KD
WB Positive control:	Hela, HepG2, MCF7 cell lysates
WB Recommended dilution:	200-1000

mlbio 梅联生物
Good elisakit producers