

兔抗 Ly6a 多克隆抗体

中文名称：兔抗 Ly6a 多克隆抗体

英文名称：Anti-Ly6a rabbit polyclonal antibody

别名：lymphocyte antigen 6 complex, locus A; TAP; Sca1; Sca-1; Ly-6A.2; Ly-6A/E; Ly-6E.1

相关类别：一抗

抗原：Ly6a

储存：冷冻 (-20℃)

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

The activity of natural killer (NK) cells is regulated by members of multiple receptor families that recognize class I MHC molecules, such as the killer cell inhibitory receptor/leukocyte immunoglobulin-like receptor (KIR/LIR) family and the C-type lectin superfamily. The Ly-6 alloantigens represent a family of phosphatidylinositol anchored proteins that play a role in the process of T lymphocyte activation. Ly-6 expression is commonly induced on T and B lymphocytes after activation by mitogens. Ly-6A/E is a phosphatidylinositol-linked transmembrane member of the Ly-6 family.

	ily that is produced in humans after 4 weeks of age. γ -interferon and TNF regulate Ly-6A/E at the level of transcription and mRNA stabilization in certain thymocytes and T-cells.
Applications:	ELISA, IHC
Name of antibody:	Ly6a
Immunogen:	Synthetic peptide of mouse Ly6a
Full name:	lymphocyte antigen 6 complex, locus A
Synonyms:	TAP; Sca1; Sca-1; Ly-6A.2; Ly-6A/E; Ly-6E.1
SwissProt:	P05533
ELISA Recommended dilution:	5000-10000
IHC positive control:	Mouse liver
IHC Recommend dilution:	20-100

