

## 兔抗 TBXAS1 多克隆抗体

中文名称：兔抗 TBXAS1 多克隆抗体

英文名称： Anti-TBXAS1 rabbit polyclonal antibody

别名： TS; TXS; CYP5; THAS; TXAS; CYP5A1; GHOSAL; BDPLT14

相关类别： 一抗

抗原： TBXAS1

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. However, this protein is considered a member of the cytochrome P450 superfamily on the basis of sequence similarity rather than functional similarity. This endoplasmic reticulum membrane protein catalyzes the conversion of prostaglandin H<sub>2</sub> to thromboxane A<sub>2</sub>, a potent vasoconstrictor and inducer of platelet aggregation. The enzyme plays a role in several pathophysiological pro

	cesses including hemostasis, cardiovascular disease, and stroke. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	TBXAS1
<b>Immunogen:</b>	Fusion protein of human TBXAS1
<b>Full name:</b>	thromboxane A synthase 1 (platelet)
<b>Synonyms :</b>	TS; TXS; CYP5; THAS; TXAS; CYP5A1; GHOSAL; BDPLT14
<b>SwissProt:</b>	P24557
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer and human breast cancer
<b>IHC Recommend dilution:</b>	50-200

