

小鼠抗 FYN 单克隆抗体

中文名称： 小鼠抗 FYN 单克隆抗体

英文名称： Anti-FYN mouse monoclonal antibody

别名： FYN proto-oncogene, Src family tyrosine kinase; SLK; SYN; p59-FYN

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： FYN

反应种属： Human, Monkey, Rat, Mouse

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:	This gene is a member of the protein-tyrosine kinase oncogene family. It encodes a membrane-associated tyrosine kinase that has been implicated in the control of cell growth. The protein associates with the p85 subunit of phosphatidylinositol 3-kinase and interacts with the fyn-binding protein. Alternatively spliced transcript variants encoding distinct isoforms exist.
Applications:	WB, IF
Name of antibody:	FYN
Immunogen:	Fusion protein of human FYN
Full name:	FYN proto-oncogene, Src family tyrosine kinase

Synonyms:	SLK; SYN; p59-FYN
SwissProt:	P06241
WB Predicted band size:	61 kDa
WB Positive control:	CEM, Ramos, Jurkat, C6, 3T3, HeLa and COS7 cell lysates
WB Recommended dilution:	500-2000
IF positive control:	HeLa cells
IF Recommend dilution:	10-50



