

## 兔抗 LY6E 多克隆抗体

中文名称：兔抗 LY6E 多克隆抗体

英文名称：Anti-LY6E rabbit polyclonal antibody

别名：RIGE; SCA2; RIG-E; SCA-2; TSA-1

抗原：LY6E

储存：冷冻（-20℃）避光

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	Induced by retinoic acid; in promyelocytic leukemia NB4 and in myeloblast HL-60 cell lines. Activated by IFN-alpha in monocytic cell line U-937 and in peripheral blood monocyte cells.
<b>Applications:</b>	IF
<b>Name of antibody:</b>	LY6E
<b>Immunogen:</b>	Synthesized peptide derived from internal of human LY6E.
<b>Full name:</b>	lymphocyte antigen 6 complex, locus E
<b>Synonyms :</b>	RIGE; SCA2; RIG-E; SCA-2; TSA-1
<b>SwissProt:</b>	Q16553

<b>IF positive control:</b>	A549 cells
<b>IF Recommend dilution:</b>	100-500

