

小鼠抗 GABPB1 单克隆抗体

中文名称： 小鼠抗 GABPB1 单克隆抗体

英文名称： Anti-GABPB1 mouse monoclonal antibody

别名： GA binding protein transcription factor beta subunit 1; E4TF1; GABPB; BABPB2; E4TF1B; GABPB2; NRF2B1; NRF2B2; GABPB-1; E4TF1-47; E4T

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： GABPB1

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:

This gene encodes the GA-binding protein transcription factor, beta subunit. This protein forms a tetrameric complex with the alpha subunit, and stimulates transcription of target genes. The encoded protein may be involved in activation of cytochrome oxidase expression and nuclear control of mitochondrial function. The crystal structure of a similar protein in mouse has been resolved as a ternary protein complex. Multiple transcript variants encoding distinct isoforms have been identified for this gene.

Applications:	WB
Name of antibody:	GABPB1
Immunogen:	Fusion protein of human GABPB1
Full name:	GA binding protein transcription factor beta subunit 1
Synonyms:	E4TF1; GABPB; BABPB2; E4TF1B; GABPB2; NRF2B1; NRF2B2; GABPB-1; E4TF1-47; E4TF1-53
SwissProt:	Q06547
WB Predicted band size:	42 KD
WB Positive control:	293T, A431, Hela, HepG2, Raji cell lysates
WB Recommended dilution:	500-2000