

## 兔抗 LY75 多克隆抗体

中文名称：兔抗 LY75 多克隆抗体

英文名称：Anti-LY75 rabbit polyclonal antibody

别名：lymphocyte antigen 75; CD205; LY-75; CLEC13B; DEC-205; GP200-MR6

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：LY75

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Acts as an endocytic receptor to direct captured antigens from the extracellular space to a specialized antigen-processing compartment (By similarity). Causes reduced proliferation of B-lymphocytes.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	LY75
<b>Immunogen:</b>	Synthetic peptide of human LY75
<b>Full name:</b>	lymphocyte antigen 75
<b>Synonyms:</b>	CD205; LY-75; CLEC13B; DEC-205; GP200-MR6

<b>SwissProt:</b>	O60449
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human breast cancer
<b>IHC Recommend dilution:</b>	25-50

