

兔抗 TCEAL3 多克隆抗体

中文名称：兔抗 TCEAL3 多克隆抗体

英文名称： Anti-TCEAL3 rabbit polyclonal antibody

别名： WEX8

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： TCEAL3

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a member of the transcription elongation factor A (SII)-like (TCEAL) gene family. Members of this family contain TFA domains and may function as nuclear phosphoproteins that modulate transcription in a promoter context-dependent manner. Multiple family members are located on the X chromosome. Alternative splicing results in multiple transcript variants encoding a single isoform.
Applications:	IHC
Name of antibody:	TCEAL3
Immunogen:	Synthesized peptide derived from internal of human

	TCEAL3/5/6.
Full name:	transcription elongation factor A (SII)-like 3
Synonyms :	WEX8
SwissProt:	Q969E4/Q5H9L2/Q6IPX3
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100

