

## 兔抗 LYL1 多克隆抗体

中文名称：兔抗 LYL1 多克隆抗体

英文名称：Anti-LYL1 rabbit polyclonal antibody

别名：bHLHa18

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：LYL1

宿主：Rabbit

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	This gene represents a basic helix-loop-helix transcription factor. The encoded protein may play roles in blood vessel maturation and hematopoiesis. A translocation between this locus and the T cell receptor beta locus (GeneID 6957) on chromosome 7 has been associated with acute lymphoblastic leukemia.
<b>Applications:</b>	IHC
<b>Name of antibody:</b>	LYL1
<b>Immunogen:</b>	Synthesized peptide derived from internal of human Lyl-1.

<b>Full name:</b>	lymphoblastic leukemia associated hematopoiesis regulator 1
<b>Synonyms :</b>	bHLHa18
<b>SwissProt:</b>	P12980
<b>IHC positive control:</b>	Human lung carcinoma tissue
<b>IHC Recommend dilution:</b>	50-100

