

兔抗 LYN 多克隆抗体

中文名称：兔抗 LYN 多克隆抗体

英文名称：Anti-LYN rabbit polyclonal antibody

别名：JTK8; p53Lyn; p56Lyn

相关类别：一抗

抗原：LYN

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a tyrosine protein kinase, which maybe involved in the regulation of mast cell degranulation, and erythroid differentiation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	LYN
Immunogen:	Fusion protein of human LYN
Full name:	v-src-1 Yamaguchi sarcoma viral related oncogene homolo

	g
Synonyms :	JTK8; p53Lyn; p56Lyn
SwissProt:	P07948
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and Human liver cancer
IHC Recommend dilution:	25-100

