

兔抗 LYPD4 多克隆抗体

中文名称：兔抗 LYPD4 多克隆抗体

英文名称：Anti-LYPD4 rabbit polyclonal antibody

别名：SMR

相关类别：一抗

储存：冷冻（-20℃）

抗原：LYPD4

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The specific function of LYPD4 is not yet known.
Applications:	ELISA, WB, IHC
Name of antibody:	LYPD4
Immunogen:	Full length fusion protein
Full name:	LY6/PLAUR domain containing 4
Synonyms :	SMR
SwissProt:	Q6UWN0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer

IHC Recommend dilution:	25-100
WB Predicted band size:	27 kDa
WB Positive control:	Human testis tissue
WB Recommended dilution:	500-2000



