

小鼠抗 GATA6 单克隆抗体

中文名称： 小鼠抗 GATA6 单克隆抗体

英文名称： Anti-GATA6 mouse monoclonal antibody

别名： GATA binding protein 6

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： GATA6

反应种属： Human

标记物： Unconjugate

克隆类型： Mouse Monoclonal

包装：
100 ul 2000

技术规格

Background:

This gene is a member of a small family of zinc finger transcription factors that play an important role in the regulation of cellular differentiation and organogenesis during vertebrate development. This gene is expressed during early embryogenesis and localizes to endo- and mesodermally derived cells during later embryogenesis and thereby plays an important role in gut, lung, and heart development. M

	utations in this gene are associated with several congenital defects.
Applications:	WB
Name of antibody:	GATA6
Immunogen:	Fusion protein of human GATA6
Full name:	GATA binding protein 6
SwissProt:	Q92908
WB Predicted band size:	60 kDa
WB Positive control:	HelaNE, Hela and 293 cell lysates
WB Recommended dilution:	200-1000

