

兔抗 LYPLA2 多克隆抗体

中文名称：兔抗 LYPLA2 多克隆抗体

英文名称：Anti-LYPLA2 rabbit polyclonal antibody

别名：APT2; APT-2; DJ886K2.4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：LYPLA2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Lysophospholipases are enzymes that act on biological membranes to regulate the multifunctional lysophospholipids. There are alternatively spliced transcript variants described for this gene but the full length nature is not known yet.
Applications:	ELISA, IHC
Name of antibody:	LYPLA2
Immunogen:	Full length fusion protein
Full name:	lysophospholipase II
Synonyms :	APT2; APT-2; DJ886K2.4

SwissProt:	O95372
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	25-100



