

小鼠抗 GC 单克隆抗体

中文名称： 小鼠抗 GC 单克隆抗体

英文名称： Anti-GC mouse monoclonal antibody

别名： GC, vitamin D binding protein; DBP; GRD3; VDBG; VDBP; GcMAF; DBP/GC; Gc-MAF; HEL-S-51

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： GC

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	The protein encoded by this gene belongs to the albumin gene family. It is a multifunctional protein found in plasma, ascitic fluid, cerebrospinal fluid and on the surface of many cell types. It binds to vitamin D and its plasma metabolites and transports them to target tissues. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	WB, IHC
Name of antibody:	GC
Immunogen:	Fusion protein of human GC

Full name:	GC, vitamin D binding protein
Synonyms:	DBP; GRD3; VDBG; VDBP; GcMAF; DBP/GC; Gc-MAF; HEL-S-51
SwissProt:	P02774
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-150
WB Predicted band size:	53 KD
WB Positive control:	Hela and T47D cell lysates
WB Recommended dilution:	500-2000