

## 小鼠抗 GIPC1 单克隆抗体

中文名称：小鼠抗 GIPC1 单克隆抗体

英文名称：Anti-GIPC1 mouse monoclonal antibody

别名：GIPC PDZ domain containing family member 1; NIP; GIPC; IIP-1; TIP-2; SEMCAP; C19orf3; Hs.6454; GLUT1CBP; RGS19IP1; SYNECTIN; SYNECTI

相关类别：一抗

储存：冷冻 (-20℃)

宿主：Mouse

抗原：GIPC1

反应种属：Human

标记物：Unconjugate

克隆类型：mouse monoclonal

### 技术规格

<b>Background:</b>	GIPC1 is a scaffolding protein that regulates cell surface receptor expression and trafficking (Lee et al., 2008 [PubMed 18775991]).[supplied by OMIM, Apr 2009]
<b>Applications:</b>	WB
<b>Name of antibody:</b>	GIPC1
<b>Immunogen:</b>	Fusion protein of human GIPC1
<b>Full name:</b>	GIPC PDZ domain containing family member 1
<b>Synonyms:</b>	NIP; GIPC; IIP-1; TIP-2; SEMCAP; C19orf3; Hs.6454; GLUT1CBP; RGS19IP1; SYNECTIN; SYNECTIN
<b>SwissProt:</b>	O14908



<b>WB Predicted band size:</b>	36 KD
<b>WB Positive control:</b>	Hela cell lysate
<b>WB Recommended dilution:</b>	500-2000