

## 兔抗 LYVE1 多克隆抗体

- 中文名称：兔抗 LYVE1 多克隆抗体
- 英文名称：Anti-LYVE1 rabbit polyclonal antibody
- 别名：HAR, XLKD1, LYVE-1, CRSBP-1
- 抗原：LYVE1
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 相关类别：一抗
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a type I integral membrane glycoprotein. The encoded protein acts as a receptor and binds to both soluble and immobilized hyaluronan. This protein may function in lymphatic hyaluronan transport and have a role in tumor metastasis.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	LYVE1
<b>Immunogen:</b>	Fusion protein of human LYVE1
<b>Full name:</b>	lymphatic vessel endothelial hyaluronan receptor 1
<b>Synonyms :</b>	HAR, XLKD1, LYVE-1, CRSBP-1

<b>SwissProt:</b>	Q9Y5Y7
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human ovarian cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	35 kDa
<b>WB Positive control:</b>	A431 cells
<b>WB Recommended dilution:</b>	100-500



