

## 小鼠抗 GLB1 单克隆抗体

中文名称: 小鼠抗 GLB1 单克隆抗体

英文名称: Anti-GLB1 mouse monoclonal antibody

别名: EBP; ELNR1; MPS4B

储存: 冷冻 (-20°C) 避光

抗原: GLB1

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes beta-galactosidase-1, a lysosomal enzyme that hydrolyzes the terminal beta-galactose from ganglioside substrates and other glycoconjugates. Defects in this gene are the cause of GM1-gangliosidosis and Morquio B syndrome. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq].
<b>Applications:</b>	WB, IF, FC
<b>Name of antibody:</b>	GLB1
<b>Immunogen:</b>	Fusion protein of human GLB1
<b>Full name:</b>	galactosidase, beta 1 (GLB1), transcript variant 2

<b>Synonyms:</b>	EBP; ELNR1; MPS4B
<b>SwissProt:</b>	P16278
<b>WB Predicted band size:</b>	76 kDa
<b>WB Positive control:</b>	Hela, SVT2, A549, COS7, MDCK, MCF-7 cell lysates
<b>WB Recommended dilution:</b>	500-2000