

小鼠抗 GLS2 单克隆抗体

中文名称： 小鼠抗 GLS2 单克隆抗体

英文名称： Anti-GLS2 mouse monoclonal antibody

别名： glutaminase 2; GA; GLS; LGA; hLGA

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： GLS2

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	The protein encoded by this gene is a mitochondrial phosphate-activated glutaminase that catalyzes the hydrolysis of glutamine to stoichiometric amounts of glutamate and ammonia. Originally thought to be liver-specific, this protein has been found in other tissues as well. Alternative splicing results in multiple transcript variants that encode different isoforms.
Applications:	WB, IHC
Name of antibody:	GLS2
Immunogen:	Fusion protein of human GLS2
Full name:	glutaminase 2

Synonyms:	GA; GLS; LGA; hLGA
SwissProt:	Q9UI32
IHC positive control:	Human liver cancer and Human kidney cancer; Human pancreas tissue
IHC Recommend dilution:	500-2000
WB Predicted band size:	66 KD
WB Positive control:	HepG2 and 293T cell lysates
WB Recommended dilution:	500-2000