

## 兔抗 M1AP 多克隆抗体

- 中文名称： 兔抗 M1AP 多克隆抗体
- 英文名称： Anti-M1AP rabbit polyclonal antibody
- 别名： meiosis 1 associated protein; D6Mm5e; C2orf65; SPATA37
- 相关类别： 一抗
- 储存： 冷冻（-20℃）
- 抗原： M1AP
- 宿主： Rabbit
- 反应种属： Human
- 标记物： Unconjugate
- 克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein that is likely to function in progression of meiosis. A similar protein in mouse plays a role in gametogenesis in both sexes. Alternate splicing results in multiple transcript variants encoding different isoforms.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	M1AP
<b>Immunogen:</b>	Fusion protein of human M1AP
<b>Full name:</b>	meiosis 1 associated protein
<b>Synonyms:</b>	D6Mm5e; C2orf65; SPATA37

<b>SwissProt:</b>	Q8TC57
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	59 kDa
<b>WB Positive control:</b>	Jurkat cell lysate
<b>WB Recommended dilution:</b>	200-1000



