

小鼠抗 GRIPAP1 单克隆抗体

中文名称： 小鼠抗 GRIPAP1 单克隆抗体

英文名称： Anti-GRIPAP1 mouse monoclonal antibody

别名： GRASP-1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Mouse

抗原： GRIPAP1

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	This gene encodes a guanine nucleotide exchange factor for the Ras family of small G proteins (RasGEF). In brain studies, the encoded protein was found with the GRIP/AMPA receptor complex. Multiple alternatively spliced transcript variants have been described that encode different protein isoforms; however, the full-length nature and biological validity of all of these variants have not been determined
Applications:	WB, IF, FC
Name of antibody:	GRIPAP1
Immunogen:	Fusion protein of human GRIPAP1
Full name:	GRIP1 associated protein 1 (GRIPAP1)

Synonyms:	GRASP-1
SwissProt:	Q4V328
WB Predicted band size:	96 kDa
WB Positive control:	Hela, A549, COS7, Jurkat, MCF-7 cell lysates
WB Recommended dilution:	500-2000