

小鼠抗 GSK3B 单克隆抗体

- 中文名称： 小鼠抗 GSK3B 单克隆抗体
- 英文名称： Anti-GSK3B mouse monoclonal antibody
- 别名： glycogen synthase kinase 3 beta
- 相关类别： 一抗
- 储存： 冷冻（-20℃）
- 宿主： Mouse
- 抗原： GSK3B
- 反应种属： Human, Mouse
- 标记物： Unconjugate
- 克隆类型： mouse monoclonal

技术规格

Background:	The protein encoded by this gene is a serine-threonine kinase belonging to the glycogen synthase kinase subfamily. It is a negative regulator of glucose homeostasis and is involved in energy metabolism, inflammation, ER-stress, mitochondrial dysfunction, and apoptotic pathways. Defects in this gene have been associated with Parkinson disease and Alzheimer disease.
Applications:	WB, IHC
Name of antibody:	GSK3B
Immunogen:	Fusion protein of human GSK3B
Full name:	glycogen synthase kinase 3 beta

SwissProt:	P49841
IHC positive control:	Human tonsil and Human colorectal cancer; Human esophagus cancer
IHC Recommend dilution:	100-500
WB Predicted band size:	47 kDa
WB Positive control:	293T, Hela cell, Mouse brain tissue lysates
WB Recommended dilution:	1000-5000





