

小鼠抗 GSTO1 单克隆抗体

中文名称： 小鼠抗 GSTO1 单克隆抗体

英文名称： Anti-GSTO1 mouse monoclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Mouse

抗 原： GSTO1

反应种属： Human

标 记 物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	The protein encoded by this gene is an omega class glutathione S-transferase (GST) with glutathione-dependent thiol transferase and dehydroascorbate reductase activities. GSTs are involved in the metabolism of xenobiotics and carcinogens. The encoded protein acts as a homodimer and is found in the cytoplasm. Three transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB
Name of antibody:	GSTO1
Immunogen:	Fusion protein of human GSTO1
Full name:	glutathione S-transferase omega 1
Synonyms:	P28; SPG-R; GSTO 1-1; GSTTLp28; HEL-S-21

SwissProt:	P78417
WB Predicted band size:	28 kDa
WB Positive control:	HepG2 whole cell lysate
WB Recommended dilution:	500-2000

