

## 小鼠抗 NOS3 单克隆抗体

中文名称： 小鼠抗 NOS3 单克隆抗体

英文名称： Anti-NOS3 mouse monoclonal antibody

别名： nitric oxide synthase 3; eNOS; ECNOS

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： NOS3

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	Nitric oxide is a reactive free radical which acts as a biologic mediator in several processes, including neurotransmission and antimicrobial and antitumoral activities. Nitric oxide is synthesized from L-arginine by nitric oxide synthases. Variations in this gene are associated with susceptibility to coronary spasm. Alternative splicing and the use of alternative promoters results in multiple transcript variants.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	NOS3
<b>Immunogen:</b>	Fusion protein of human NOS3
<b>Full name:</b>	nitric oxide synthase 3
<b>Synonyms:</b>	eNOS; ECNOS

<b>SwissProt:</b>	P29474
<b>WB Predicted band size:</b>	133 kDa
<b>WB Positive control:</b>	Human heart tissue lysate
<b>WB Recommended dilution:</b>	1000-5000

