

兔抗 TEKT3 多克隆抗体

中文名称：兔抗 TEKT3 多克隆抗体

英文名称：Anti-TEKT3 rabbit polyclonal antibody

抗 原： TEKT3

储 存： 冷冻 (-20℃)

宿 主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene product belongs to the tektin family of proteins. Tektins comprise a family of filament-forming proteins that are coassembled with tubulins to form ciliary and flagellar microtubules. The exact function of this gene is not known. Structural component of ciliary and flagellar microtubules. Forms filamentous polymers in the walls of ciliary and flagellar microtubules. Required for progressive sperm mobility.
Applications:	ELISA, IHC
Name of antibody:	TEKT3
Immunogen:	Fusion protein of human TEKT3
Full name:	Tektin 3
SwissProt:	Q9BXF9

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human liver cancer
IHC Recommend dilution:	25-100

