

## 小鼠抗 GST 单克隆抗体

中文名称： 小鼠抗 GST 单克隆抗体

英文名称： Anti-GST mouse monoclonal antibody

相关类别： 标签抗体

储 存： 冷冻（-20℃） 避光

宿 主： Mouse

抗 原： GST

反应种属： GST tag

标 记 物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	GST (Glutathione S-Transferase) is a 26kDa protein which was widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins. GST isoenzymes appear to play a central role in the parasite detoxification system. Other functions are also suspected including a role in increasing the solubility of haematin in the parasite gut.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	GST
<b>Immunogen:</b>	Recombinant GST protein
<b>Full name:</b>	Glutathione S-Transferase
<b>WB Positive control:</b>	six different HEK293T over-expression lysates containing different GST tagged recombinant proteins, and one GST-D



	DK bacterial lysate produced in E.coli.
<b>WB Recommended dilution:</b>	2000-10000